

Medical Examination Nitrile Gloves

KS-ST RT021(Multi-Language)

Our gloves are manufactured by overcoming numerous tests based on the corresponding regulations. KINGFA has a daily production capacity of 1,000,000 units per line. At the end of 2021, KINGFA will have 96 functional production lines, reaching a total daily production capacity of over 100,000,000 units.



PRODUCT INFORMATION

Material	100% Nitrile		
Color	Blue		
Cuff length	Standard		
Glove length(mm/inches) min	240 / 9.5		
Powder content	Powder-Free		
External glove surface	Textured Finger		
Freedom from holes(Inspection level I)	1.5 AQL		
Palm thickness (mm/mil)	≥0.05 / 2.0		
Finger thickness (mm/mil)	≥0.08 / 3.1		
Physical Requirement (median value)	Before Aging	Breakage Force(N)	≥6N
	After Aging	Breakage Force(N)	≥6N
Application Temperature	≤40°C		
Size	Median Length(mm)	Median Width(mm)	Unit Weight(gram)
S	≥240	80±10	3.5±0.3
M	≥240	95±10	3.7±0.3
L	≥240	110±10	4.0±0.3
XL	≥240	≥110	4.3±0.3

INTENDED USE

Industries:

- Pharmaceutical
- Dentistry and hygiene
- Healthcare
- Pet care
- Food processing

FEATURES

- Excellent mechanical strength provides a high level of hand protection.
- Textured in finger tips for a secure grip.
- Protection against bacteria and fungi.
- Skin irritation&sensitization tested.
- Multi-Language Packaging: EN,DE,FR,NL,IT,ES.

STANDARDS AND CERTIFICATIONS

EN 455,EN 374,EN 1186
EU Module B+C2 Certificate

MANUFACTURING ACCREDITATIONS

ISO/IEC 17025, ISO 9001, ISO 13485
OHSAS 18001, ISO 14001. IECQ QC 080000

REGULATORY COMPLIANCE

MDD 93/42/EEC,REACH,EU 2016/425,
EC 10/2011,EC 1935/2004

PACKAGING

100 pcs/box
Box Size:225*120*63 mm
10 boxes/carton
Carton Size:330*250*240 mm

CONTAINER LOADING

(For Reference Only/Without Pallets)

20GP:1500 cartons
40GP:3100 cartons
40HQ:3464 cartons

